

RRAGA Polyclonal Antibody



Catalog Number: 165527

Product name

RRAGA Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to RRAGA

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-313 of human RRAGA (NP_006561.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

Validations



Immunofluorescence - RRAGA Polyclonal Antibody

Immunofluorescence analysis of A549 cells using RRAGA antibody . Blue: DAPI for nuclear staining.