

ADRA1B Polyclonal Antibody



Catalog Number: 165578

Product name

ADRA1B Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ADRA1B

Preparation

Antigen: A synthetic peptide of human ADRA1B

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations



Western blot - ADRA1B Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ADRA1B antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Immunofluorescence - ADRA1B Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using ADRA1B antibody at dilution of 1:100. Blue: DAPI for nuclear staining.