

STX1A Polyclonal Antibody



Catalog Number: 165618

Product name

STX1A Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to STX1A

Preparation

Antigen: A synthetic peptide of human STX1A

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:100 - 1:200


IF 1:50 - 1:200

Validations




Immunohistochemistry - STX1A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using STX1A antibody (40x lens).



Immunohistochemistry - STX1A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using STX1A antibody (40x lens).



Immunohistochemistry - STX1A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded

STX1A Polyclonal Antibody



Catalog Number: 165618

mouse brain using STX1A antibody (40x lens).