

# JMJD4 Polyclonal Antibody



Catalog Number: 165652

## Product name

Rabbit IgG

JMJD4 Polyclonal Antibody

## Applications

## Specificity

WB

Human, Mouse, Rat

## Dilutions

## Antibody description


WB 1:1000 - 1:5000

Polyclonal antibody to JMJD4

## Validations

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-463 of human JMJD4 (NP\_075383.2).



## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Western blot - JMJD4 Polyclonal Antibody

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

Western blot analysis of extracts of various cell lines, using JMJD4 antibody at 1:5000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

## Clonality

Polyclonal

## Ig Type