

RFC5 Polyclonal Antibody



Catalog Number: 165676

General Information

Product name:

RFC5 Polyclonal Antibody

Immunogen:

Applications:

WB, IHC

Specificity:

Human, Mouse, Rat

Antibody description:

Polyclonal antibody to RFC5

Properties

Preparation:

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 181-340 of human RFC5 (NP_031396.1).

Formulation:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage:

Store at -20°C. Avoid freeze / thaw cycles.

Clonality:

Polyclonal


Ig Type:Rabbit IgG

Applications

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - RFC5 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using RFC5 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Immunohistochemistry - RFC5 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using RFC5 antibody at dilution of 1:100 (40x lens).