

# Gm13125 Polyclonal Antibody



Catalog Number: 165699

## Product name

Gm13125 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to Gm13125

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 67-274 of mouse Gm13125 (NP\_001108549.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:100 - 1:200

## Validations



## Western blot - Gm13125 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using Gm13125 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



## Immunohistochemistry - Gm13125 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human liver cancer using Gm13125 antibody at dilution of 1:100 (40x lens).

# Gm13125 Polyclonal Antibody



Catalog Number: 165699

---



Immunohistochemistry - Gm13125 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using Gm13125 antibody at dilution of 1:100 (40x lens).