

DAPK2 Polyclonal Antibody



Catalog Number: 165743

Product name

DAPK2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to DAPK2

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 300 to the C-terminus of human DAPK2 (NP_055141.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

Validations



Immunofluorescence - DAPK2 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using DAPK2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.