

UBE2U Polyclonal Antibody



Catalog Number: 165831

Product name

UBE2U Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to UBE2U

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-226 of human UBE2U (NP_689702.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:100

Validations




Western blot - UBE2U Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using UBE2U antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Immunohistochemistry - UBE2U Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver using UBE2U antibody at dilution of



1:100 (40x lens).

UBE2U Polyclonal Antibody



Catalog Number: 165831

Immunohistochemistry - UBE2U Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using UBE2U antibody at dilution of 1:100 (40x lens).