

# ALAD Polyclonal Antibody



Catalog Number: 165844

## Product name

ALAD Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to ALAD

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human ALAD (NP\_000022.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations



Western blot - ALAD Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ALAD antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - ALAD Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using ALAD antibody at dilution of 1:100 (40x lens).