

# GNRH2 Polyclonal Antibody



Catalog Number: 165861

## Product name

GNRH2 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to GNRH2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 24-120 of human GNRH2 (NP\_001492.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

## Validations




Western blot - GNRH2 Polyclonal Antibody

Western blot analysis of extracts of 293T cells, using GNRH2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - GNRH2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using GNRH2 antibody at dilution of 1:100



(40x lens).

# GNRH2 Polyclonal Antibody




Catalog Number: 165861

---


Immunohistochemistry - GNRH2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human esophagus using GNRH2 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - GNRH2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using GNRH2 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - GNRH2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using GNRH2 antibody at dilution of 1:100 (40x lens).