

# NDUFA7 Polyclonal Antibody



Catalog Number: 165871

## Product name

NDUFA7 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to NDUFA7

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human NDUFA7 (NP\_004992.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:100

## Validations



Western blot - NDUFA7 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using NDUFA7 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 30s.

Western blot - NDUFA7 Polyclonal Antibody


Western blot analysis of extracts from normal (control) and NDUFA7 knockout 293T cells, using NDUFA7 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time: 5s.

# NDUFA7 Polyclonal Antibody



Catalog Number: 165871

---



Immunohistochemistry - NDUFA7 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver injury using NDUFA7 antibody at dilution of 1:100 (40x lens).