

INADL Polyclonal Antibody



Catalog Number: 165898

Product name

INADL Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to INADL

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 450-680 of human INADL (NP_795352.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:100

Validations




Western blot - INADL Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using INADL antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.



Immunohistochemistry - INADL Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using INADL antibody at dilution of



1:100 (40x lens).


INADL Polyclonal Antibody



Catalog Number: 165898

Immunohistochemistry - INADL Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver injury using INADL antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - INADL Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using INADL antibody at dilution of 1:100 (40x lens).