

# ARID2 Polyclonal Antibody

Catalog Number: 165979

## Product name

ARID2 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to ARID2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1586-1835 of human ARID2 (NP\_689854.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions


WB 1:500 - 1:2000

## Validations



Western blot - ARID2 Polyclonal Antibody

Western blot analysis of extracts of LO2 cells, using ARID2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 30s.

Western blot - ARID2 Polyclonal Antibody

Western blot analysis of extracts from normal (control) and ARID2 knockout HeLa cells, using ARID2 antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time: 10s.