

SOX10 Polyclonal Antibody



Catalog Number: 165990

General Information

Product name:

SOX10 Polyclonal Antibody

Immunogen:

Applications:

WB, IHC

Specificity:

Human, Mouse, Rat

Antibody description:

Polyclonal antibody to SOX10

Properties

Preparation:

Antigen: A synthetic Peptide of human SOX10

Formulation:

PBS with 0.02% sodium azide, pH7.3.

Storage:

Store at 4°C. Avoid freeze / thaw cycles.

Clonality:

Polyclonal


Ig Type:Rabbit IgG

Applications

WB 1:500 - 1:1000

IHC 1:20 - 1:50

Validations



Western blot - SOX10 Polyclonal Antibody

Western blot analysis of extracts of A-375 cells, using SOX10 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.