

GSTM1 Polyclonal Antibody



Catalog Number: 166019

Product name

GSTM1 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to GSTM1

Preparation

Antigen: A synthetic peptide of human GSTM1

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, FC


Dilutions

WB 1:200 - 1:500

IHC 1:50 - 1:100

FC 1:20 - 1:50

Validations



Western blot - GSTM1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using GSTM1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.