

# SIGLEC6 Polyclonal Antibody



Catalog Number: 166037

## Product name

SIGLEC6 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to SIGLEC6

## Preparation

Antigen: A synthetic Peptide of human SIGLEC6

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:1000

IHC 1:20 - 1:50

## Validations



Western blot - SIGLEC6 Polyclonal Antibody

Western blot analysis of extracts of mouse spleen, using SIGLEC6 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.