

# CCDC6 Polyclonal Antibody



Catalog Number: 166046

## Product name

CCDC6 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to CCDC6

## Preparation

Antigen: A synthetic Peptide of human CCDC6

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:1000

## Validations



Western blot - CCDC6 Polyclonal Antibody

Western blot analysis of extracts of K-562 cells, using CCDC6 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.