

FGF5 Polyclonal Antibody



Catalog Number: 166089

Product name

Rabbit IgG

FGF5 Polyclonal Antibody

Applications

Specificity

WB

Human

Dilutions

Antibody description


WB 1:200 - 1:2000

Polyclonal antibody to FGF5

Validations

Preparation

Antigen: A synthetic peptide of human FGF5



Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Western blot - FGF5 Polyclonal Antibody

Clonality

Polyclonal

Western blot analysis of extracts of A-375 cells, using FGF5 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.

Ig Type