

# SLC24A4 Polyclonal Antibody



Catalog Number: 166096

## Product name

SLC24A4 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to SLC24A4

## Preparation

Antigen: A synthetic peptide of human SLC24A4

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB

## Dilutions

WB 1:500 - 1:1000

## Validations

Human lung



Western blot - SLC24A4 Polyclonal Antibody

Western blot analysis of extracts of human lung, using SLC24A4 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.