

DAB2IP Monoclonal Antibody



Catalog Number: 166134

Product name

DAB2IP Monoclonal Antibody

Specificity

Human, Mouse

Antibody description

Monoclonal antibody to DAB2IP

Preparation

Antigen: A synthetic peptide of human DAB2IP

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG

Applications


WB, IP

Dilutions

WB 1:500 - 1:2000

IP 1:50 - 1:100

Validations



Western blot - DAB2IP Monoclonal Antibody

Western blot analysis of extracts of HeLa cells, using DAB2IP antibody. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.