

MT-CYB Polyclonal Antibody



Catalog Number: 166155

Product name

MT-CYB Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to MT-CYB

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse MT-CYB (NP_904340.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions


WB 1:200 - 1:1000

Validations



Western blot - MT-CYB Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using MT-CYB antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 1s.

Western blot - MT-CYB Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using MT-CYB antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit .Exposure time: 5s.