

COL9A3 Polyclonal Antibody



Catalog Number: 166159

Product name

COL9A3 Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to COL9A3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 445-684 of human COL9A3 (NP_001844.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions

IHC 1:100 - 1:200

Validations



Immunohistochemistry - COL9A3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse heart using COL9A3 antibody at dilution of 1:100 (40x lens).