

# HOXA13 Polyclonal Antibody



Catalog Number: 166174

## Product name

HOXA13 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to HOXA13

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 140-320 of human HOXA13 (NP\_000513.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:1000 - 1:4000

IF 1:50 - 1:200

## Validations



Western blot - HOXA13 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HOXA13 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.