

FYB Monoclonal Antibody



Catalog Number: 166230

Product name

FYB Monoclonal Antibody

Specificity

Human

Antibody description

Monoclonal antibody to FYB

Preparation

Antigen: Recombinant protein of human FYB

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG


Applications

WB

Dilutions

WB 1:200 - 1:1000

Validations



Western blot - FYB Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using FYB antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.