

TRIM33 Monoclonal Antibody



Catalog Number: 166235

Product name

TRIM33 Monoclonal Antibody

Specificity

Human

Antibody description

Monoclonal antibody to TRIM33

Preparation

Antigen: Recombinant protein of human TRIM33

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - TRIM33 Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using TRIM33 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.