

TMED1 Polyclonal Antibody



Catalog Number: 166260

Product name

Rabbit IgG

TMED1 Polyclonal Antibody

Applications

Specificity

WB

Human, Mouse

Dilutions

Antibody description

WB 1:200 - 1:500

Polyclonal antibody to TMED1

Validations

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 24-194 of human TMED1 (NP_006849.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type



Western blot - TMED1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TMED1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.