

# Mouse anti GST-Tag mAb



Catalog Number: 166288

## Product name

Mouse anti GST-Tag mAb

## Antibody description

Monoclonal antibody to Mouse anti GST-Tag

## Preparation

Antigen: Recombinant protein of GST tag

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG

## Applications

WB, IF, IP


## Dilutions

WB 1:2000 - 1:5000

IF 1:50 - 1:200

IP 1:20 - 1:50

## Validations



Western blot - Mouse anti GST-Tag mAb

Western blot analysis of over-expressed GST fusion protein in 293 cell using GST-Tag antibody at different dilution. Each lane was loaded with 2 ug cell lysate. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.