

# Mouse anti Myc-Tag mAb



Catalog Number: 166292

## Product name

Mouse anti Myc-Tag mAb

## Antibody description

Monoclonal antibody to Mouse anti Myc-Tag

## Preparation

Antigen: A synthetic peptide of Myc-Tag

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG

## Applications

WB, IF, IP


## Dilutions

WB 1:2000 - 1:5000

IF 1:50 - 1:200


IP 1:50 - 1:100

## Validations



Western blot - Mouse anti Myc-Tag mAb

Western blot analysis of over-expressed Myc-tagged protein in 293T cells using Mouse anti Myc-Tag antibody at different dilution. Each lane was loaded with 2 ug cell lysate. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.




Western blot - Mouse anti Myc-Tag mAb

Western blot analysis of over-expressed Myc-tagged protein in 293T cells using Mouse anti Myc-Tag antibody. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

# Mouse anti Myc-Tag mAb


Catalog Number: 166292



Basic Kit .Exposure time: 30s.

Immunofluorescence - Mouse anti Myc-Tag mAb

Immunofluorescence analysis of 293T cells using Mouse anti Myc-Tag mAb . Blue: DAPI for nuclear



staining.

Immunoprecipitation - Mouse anti Myc-Tag mAb

Immunoprecipitation of over-expressed Myc-tagged protein in 293T cells incubated using Myc-tag antibody . A mock served as negative control and over-expressed 293T cell lysate served as positive control.