

# Rabbit anti GFP-Tag pAb



Catalog Number: 166296

## Product name

Rabbit anti GFP-Tag pAb

## Antibody description

Polyclonal antibody to Rabbit anti GFP-Tag

## Preparation

Antigen: Recombinant protein of GFP

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:2000 - 1:3000

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations



Western blot - Rabbit anti GFP-Tag pAb

Western blot analysis of extracts of normal 293T cells and 293T transfected with EGFP, using GFP antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.



Immunofluorescence - Rabbit anti GFP-Tag pAb

Immunofluorescence analysis of GFP transgenic 293T cells using GFP-Tag antibody. Green: GFP expression. Blue: DAPI for nuclear staining.