

Mouse anti TAP-Tag mAb



Catalog Number: 166298

Product name

Mouse anti TAP-Tag mAb

Antibody description

Monoclonal antibody to Mouse anti TAP-Tag

Preparation

Antigen:

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG

Applications


WB, IP

Dilutions

WB 1:5000 - 1:10000

IP 1:100 - 1:200

Validations



Western blot - Mouse anti TAP-Tag mAb

Western blot analysis of recombinant TAP protein with mouse anti TAP-Tag mAb at dilution of 1:5000. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.