

# Mouse anti MBP-Tag mAb



Catalog Number: 166300

## Product name

Mouse anti MBP-Tag mAb

## Antibody description

Monoclonal antibody to Mouse anti MBP-Tag

## Preparation

Antigen: Recombinant protein of MBP-Tag

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG


## Applications

WB

## Dilutions

WB 1:5000 - 1:10000

## Validations



Western blot - Mouse anti MBP-Tag mAb

Western blot analysis of recombinant MBP protein with Mouse anti MBP-Tag mAb at dilution of 1:5000.