

HRP-conjugated Mouse anti Myc-Tag mAb



Catalog Number: 166303

Product name

HRP-conjugated Mouse anti Myc-Tag mAb

Applications

WB

Antibody description

Monoclonal antibody to HRP-conjugated Mouse anti Myc-Tag

Dilutions

WB 1:2000 - 1:5000

Preparation

Antigen: A synthetic Peptide of human Myc-Tag

Formulation

PBS, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.


Clonality

Monoclonal

Ig Type

Mouse

Validations



Western blot - HRP-conjugated Mouse anti Myc-Tag mAb

Western blot analysis of recombinant Myc protein with HRP-conjugated Mouse anti Myc-Tag mAb at dilution of 1:5000. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.