

# Rabbit anti V5-Tag pAb



Catalog Number: 166305

## Product name

Rabbit anti V5-Tag pAb

## Antibody description

Polyclonal antibody to Rabbit anti V5-Tag

## Preparation

Antigen: A synthetic peptide of V5 epitope

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:1000 - 1:2000

## Validations



Western blot - Rabbit anti V5-Tag pAb

Western blot analysis of over-expressed V5-tagged protein in 293T cell using V5-tag antibody at different dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.