

HRP-conjugated Mouse anti HA-Tag mAb



Catalog Number: 166306

Product name

HRP-conjugated Mouse anti HA-Tag mAb

Antibody description

Monoclonal antibody to HRP-conjugated Mouse anti HA-Tag

Preparation

Antigen: A synthetic Peptide of human HA-Tag

Formulation

PBS, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse


Applications

WB

Dilutions

WB 1:1000 - 1:3000

Validations



Western blot - HRP-conjugated Mouse anti HA-Tag mAb

Western blot analysis of recombinant HA protein with HRP-conjugated Mouse anti HA-Tag mAb at dilution of 1:2000. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.