

# Phospho-MAP3K5-S83 pAb



Catalog Number: 166321

## Product name

Phospho-MAP3K5-S83 pAb

## Specificity

Human

## Antibody description

Polyclonal antibody to Phospho-MAP3K5-S83

## Preparation

Antigen: A phospho specific peptide corresponding to residues surrounding S83 of human MAP3K5

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - Phospho-MAP3K5-S83 pAb

Western blot analysis of extracts of various cell lines, using Phospho-MAP3K5-S83 antibody at 1:1000 dilution. K562 cell lysate were treated by CIP (20ul CIP for each 400ul cell lysate) at 37u2103 for 1 hour, treated by TNF-u03b1 (20ng/ml) for 30 minutes after serum-starvation overnight or treated by ATP (5mM) at 30u2103 for 1 hour. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA. Detection: ECL Basic Kit. Exposure time: 60s.