

# Phospho-MAP3K5-T918 pAb



Catalog Number: 166323

## Product name

Phospho-MAP3K5-T918 pAb

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to Phospho-MAP3K5-T918

## Preparation

Antigen: A synthetic phosphorylated peptide around T918 of human MAP3K5 .

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - Phospho-MAP3K5-T918 pAb

Western blot analysis of extracts of HeLa cells, using Phospho-MAP3K5-T918 antibody at 1:1000 dilution. HeLa cells were treated by UV for 15-30 minutes. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA.