

# Phospho-IL10RA-S319/S323 pAb



Catalog Number: 166414

## General Information

### Product name:

Phospho-IL10RA-S319/S323 pAb

### Immunogen:

### Applications:

WB

### Specificity:

Human

### Antibody description:

Polyclonal antibody to Phospho-IL10RA-S319/S323

## Properties

### Preparation:

Antigen: A phospho specific peptide corresponding to residues surrounding S319/S323 of human IL10RA

### Formulation:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

### Storage:

Store at -20°C. Avoid freeze / thaw cycles.

### Clonality:


Polyclonal

**Ig Type:**Rabbit IgG

## Applications

WB 1:500 - 1:2000

## Validations



Western blot - Phospho-IL10RA-S319/S323 pAb

Western blot analysis of extracts from HeLa cells, using Phospho-IL10RA-S319/S323 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% BSA.