

Phospho-MARCKS-S170 pAb



Catalog Number: 166432

General Information

Product name:

Phospho-MARCKS-S170 pAb

Immunogen:

Applications:

WB, IHC

Specificity:

Human

Antibody description:

Polyclonal antibody to Phospho-MARCKS-S170

Properties

Preparation:

Antigen: A phospho specific peptide corresponding to residues surrounding S170 of human MARCKS

Formulation:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage:

Store at -20°C. Avoid freeze / thaw cycles.

Clonality:

Polyclonal


Ig Type:Rabbit IgG

Applications

WB 1:500 - 1:2000


IHC 1:50 - 1:100

Validations



Western blot - Phospho-MARCKS-S170 pAb

Western blot analysis of extracts from MCF7 cells, using phospho-MARCKS-S170 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% BSA.



Immunohistochemistry - Phospho-MARCKS-S170 pAb

Immunohistochemistry of paraffin-embedded human lung carcinoma using Phospho-MARCKS-S170 antibody .