

Phospho-PRKCZ-T410 pAb



Catalog Number: 166444

General Information

Product name:

Phospho-PRKCZ-T410 pAb

Immunogen:

Applications:

IHC

Specificity:

Human, Mouse, Rat

Antibody description:

Polyclonal antibody to Phospho-PRKCZ-T410

Properties

Preparation:

Antigen: A phospho specific peptide corresponding to residues surrounding T410 of human PRKCZ

Formulation:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage:

Store at -20°C. Avoid freeze / thaw cycles.

Clonality:


Polyclonal

Ig Type:Rabbit IgG

Applications

IHC 1:50 - 1:100

Validations



Immunohistochemistry - Phospho-PRKCZ-T410 pAb

Immunohistochemistry of paraffin-embedded human lung carcinoma tissue, using Phospho-PRKCZ-T410 antibody .