

Phospho-FYN-Y416/Lck-Y416 pAb



Catalog Number: 166465

General Information

Product name:

Phospho-FYN-Y416/Lck-Y416 pAb

Immunogen:

Applications:

WB

Specificity:

Human, Mouse, Rat

Antibody description:

Polyclonal antibody to Phospho-FYN-Y416/Lck-Y416

Properties

Preparation:

Antigen: A phospho specific peptide corresponding to residues surrounding Y416 of human FYN and Y416 of human Lck

Formulation:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage:

Store at -20°C. Avoid freeze / thaw cycles.

Clonality:


Polyclonal

Ig Type:Rabbit IgG

Applications

WB 1:500 - 1:2000

Validations



Western blot - Phospho-FYN-Y416/Lck-Y416 pAb

Western blot analysis of extracts of NIH3T3 cells, using Phospho-FYN-Y416/Lck-Y416 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% BSA.