

Phospho-VASP-S157 pAb



Catalog Number: 166510

General Information

Product name:

Phospho-VASP-S157 pAb

Immunogen:

Applications:

WB, IHC, IP

Specificity:

Human, Mouse, Rat

Antibody description:

Polyclonal antibody to Phospho-VASP-S157

Properties

Preparation:

Antigen: A synthetic phosphorylated peptide around S157 of human VASP (NP_003361.1).

Formulation:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage:

Store at -20°C. Avoid freeze / thaw cycles.

Clonality:

Polyclonal

Ig Type:Rabbit IgG


Applications

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IP 1:50 - 1:100

Validations



Western blot - Phospho-VASP-S157 pAb

Western blot analysis of extracts of A431 cells, using Phospho-VASP-S157 antibody at 1:2000 dilution. A431 cells were treated by Forskolin (10uM) for 30 minutes after serum-starvation overnight. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA. Detection: ECL Basic Kit. Exposure time: 1s.



Immunoprecipitation - Phospho-VASP-S157 pAb

Immunoprecipitation analysis of 200ug extracts of A-431 cells, using 3 ug Phospho-VASP-S157 pAb. Western blot was performed from the immunoprecipitate using Phospho-VASP-S157 pAb at a dilution of 1:1000. A-431 cells were treated by Forskolin (10 uM) at 37u2103 for 30 minutes after serum-starvation overnight.