

EIF2AK3 Monoclonal Antibody



Catalog Number: 166607

General Information

Product name:

EIF2AK3 Monoclonal Antibody

Immunogen:

Applications:

WB

Specificity:

Human, Mouse, Rat

Antibody description:

Monoclonal antibody to EIF2AK3

Properties

Preparation:

Antigen: Recombinant protein of human EIF2AK3

Formulation:

PBS with 0.02% sodium azide, pH7.3.

Storage:

Store at 4°C (regular) or -20°C (long term). Avoid freeze / thaw cycles.

Clonality:


Monoclonal

Ig Type: Mouse IgG

Applications

WB 1:500 - 1:1000

Validations



Western blot - EIF2AK3 Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using EIF2AK3 antibody. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.