

Anti-Caspase-3 Rabbit antibody



Catalog Number: 166660

Product name

Anti-Caspase-3 Rabbit antibody

Specificity

Human, Mouse

Antibody description

Caspase-3 Rabbit polyclonal antibody

Preparation

Antigen: Purified recombinant human Caspase-3 protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB


Dilutions

WB: 1/1000

Validations



Western blot detection of Caspase-3 in Jurkat and 3T3 cell lysates using Caspase-3 Rabbit pAb (1:1000 diluted). Predicted band size: 35KDa.



Observed band size:35KDa.

Western blot analysis of extracts from Jurkat cells, untreated (-) or treated with Staurosporine (1 μ M; 3 hr), using Caspase-3 Rabbit pAb (166660, upper) or β -Actin (8F10) Mouse mAb (lower).