

Anti-LC3A/B Rabbit antibody



Catalog Number: 166671

Product name

Anti-LC3A/B Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

LC3A/B Rabbit polyclonal antibody

Preparation

Antigen: Synthetic peptide.

Formulation

Purified Rabbit polyclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB


Dilutions

WB: 1/1000

Validations




Western blot detection of LC3A/B in Rat Hippocampus, Rat Cerebral Cortex and Rat Brain lysates using LC3A/B Rabbit pAb (1:1000 diluted). Predicted band size: 15KDa. Observed band



size:14, 16KDa.

Western blot detection of LC3A/B in 3T3 and Rat Brain lysates using LC3A/B Rabbit pAb (1:1000 diluted). Predicted band size: 15KDa. Observed



band size:14, 16KDa.

Western blot detection of LC3A/B in Rat Brain, C2C12, 3T3 and HeLa cell lysates using LC3A/B Rabbit pAb (1:1000 diluted). Predicted band size: 15KDa. Observed band size:14, 16KDa.