

Anti-GABARAP Rabbit antibody



Catalog Number: 166697

Product name

Anti-GABARAP Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

GABARAP Rabbit polyclonal antibody

Preparation

Antigen: Synthetic peptide.

Formulation

Purified Rabbit polyclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB


Dilutions

WB: 1/1000

Validations



Western blot detection of GABARAP in A549, 3T3 and HeLa cell lysates using GABARAP Rabbit pAb (1:1000 diluted). Predicted band size: 14KDa. Observed band size:14, 16KDa.



Western blot detection of GABARAP in Rat Hippocampus, Rat Cerebral Cortex, Rat Brain and A549 cell lysates using GABARAP Rabbit pAb (1:1000 diluted). Predicted band size: 14KDa. Observed band size:14, 16KDa.