

Anti-Bcl-2 Rabbit antibody

Catalog Number: 166699

Product name

Anti-Bcl-2 Rabbit antibody

Specificity

Human, Monkey

Antibody description

Bcl-2 Rabbit polyclonal antibody

Preparation

Antigen: Peptide sequence around aa. 54~58 (G-H-T-P-H) derived from human Bcl2.

Formulation

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB


Dilutions

WB::1:300

Validations



Western blot detection of Bcl2 in a549, HeLa and COS7 cell lysates using Bcl2 rabbit pAb (dilution 1:300). Predicted band size: 26 kDa. Observed band



size: 26 kDa.

Western blot detection of Bcl2 in HeLa, Jurkat and PC-12 cell lysates using Bcl2 rabbit pAb (1:1000 diluted). Predicted band size: 26 kDa. Observed band size: 26 kDa.