

# Anti-Cdk4 Rabbit antibody



Catalog Number: 166706

## General Information

### Product name:

Anti-Cdk4 Rabbit antibody

### Immunogen:

### Applications:

WB, IHC-P

### Specificity:

Human

### Antibody description:

Cdk4 Rabbit polyclonal antibody

## Properties

### Preparation:

Antigen: Peptide sequence around aa.295~299(L-H-K-D-E)derived from Human CDK4.

### Formulation:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

### Storage:

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

### Clonality:


Polyclonal

Ig Type:IgG


## Applications

WB:1/1000IHC:1/50-1/100

## Validations



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using CDK4 Antibody .



Western blot analysis of extracts from HeLa, 293T and HUVEC cells using CDK4 rabbit pAb (1:1000 diluted). Predicted band size: 34 kDa. Observed band size: 34 kDa.