

Anti-NSE Rabbit antibody



Catalog Number: 166724

Product name

Rabbit IgG

Anti-NSE Rabbit antibody

Applications

Specificity

WB, FC, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000-1/5000

NSE Rabbit polyclonal antibody

FC: 1/20

Preparation

IP: 1/20

Antigen: Recombinant protein of human NSE

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type



Western blot detection of NSE in Human Fetal Brain lysates using NSE Rabbit pAb(1:1000 diluted). Predicted band size: 47KDa. Observed band size: 47KDa.