

# Anti-BDNF Rabbit antibody



Catalog Number: 166749

## Product name

Anti-BDNF Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

BDNF Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human BDNF

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P, ICC/IF

## Dilutions


WB: 1/1000-1/2000

IHC: 1/50-1/200

## Validations




Immunofluorescent analysis of HeLa cells, using



BDNF Antibody.

Western blot analysis of extracts of Mouse heart lysate, using BDNF antibody



Immunohistochemical analysis of paraffin-embedded human brain, using BDNF Antibody.