

Anti-CD9 Rabbit antibody



Catalog Number: 166764

Product name

Anti-CD9 Rabbit antibody

Specificity

Human

Antibody description

CD9 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human CD9

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC, IP

Dilutions


WB: 1/1000

IHC: 1/200

FC: 1/100

IP: 1/100

Validations



Western blot detection of CD9 in HeLa cell lysates using CD9 Rabbit pAb(1:1000 diluted).Predicted band size:25kDa.Observed band size:25kDa.